

# Take care of your dentist's health first.




Dental clinics are always contaminated by sprays (aerosols) smaller than 5  $\mu\text{m}$  from dental equipment such as hand-pieces, ultrasonic calculus removers, and air water syringes that are used on a daily basis. It remains in the air for more than 24 hours and can be inhaled into the lungs and can become infected by splashing microscopic particles on the dentist's face and the patient's chest during treatment.

Since the **Clean-D** dental mobile dust collector is sucked under strong vacuum pressure directly during treatment, it is possible to prevent cross-infection caused by microorganisms generated during treatment.

## Clean D Dental Mobile Dust Collector

Model Name	Clean D
Main Power	230V, 50/60Hz
Consumption Power	1.38KW
Size (L x W x H)	37 x 34 x 120 cm
Weight	28 $\pm$ 0.5 kg
Expendables	HEPA + Carbon Filter
Auxiliary components	Extraction Arm



3 Layer Filter



Prevent suction of medical tools